## Spyder-R P326AR Rehab

## PATHMASTER

## Powerchair

Eclipse Pathmaster Series – Freedom to move where you want when you want. Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.

- Pathmaster Spyder-R 10 inch mid-wheel powerchair
- > Full suspension
- Off-board 4A charger with 36AH batteries
- > 16" 17" 18" width and depth adjustable
- > Tension adjustable back
- > Swing-away joystick
- > Colours: red or blue







## P326AR Rehab Specifications

Total wt. w/bat. Safe climbo Weight capacity **Dimensions** LXWXH Motor **Battery** Tire Driving dist. Turning radius Speed 34"x 24"x 40" 88kg/194LBS 2 х .7нр 8крн/5мрн 2 х 36ан F6" M10" R6" 28км/18м 50cm/20" 6°/10% 136kg/300lbs All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style