Spyder XL-R P327R Rehab

PATHMASTER

Powerchair

Eclipse Pathmaster Series – Freedom to move where you want when you want. Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.

- Pathmaster Spyder XL-R 12 inch mid-wheel powerchair
- > Full suspension
- > Off-board 5A charger with 55AH batteries
- > 18" 19" 20" width and depth adjustable
- > Tension adjustable back
- > Swing-away joystick
- > Colours: red or blue







P327R Rehab Specifications

Total wt. w/bat. Safe climbo Weight capacity **Dimensions** LXWXH Motor Tire Driving dist. **Turning radius** Speed **Battery** 35"x 27"x 39" 110kg/242LBS 2 х .8нр 8крн/5мрн 2 х 55ан F6" M12" R6" 32км/20м 53cm/21" 8°/14% 205kg/450lbs All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style