Bigfoot S846

TRAILMASTER

Full-size scooter

The Eclipse Trailmaster series allows individuals to experience the outdoors in feature rich design, style, performance & comfort. They also come standard with an LCD digital dashboard. This series is for those individuals who want to travel a little further, in style. The full-size series are the most popular scooters in our portfolio.



> Digital LCD dashboard

> Full light package

> Front & rear suspension w/ 400 lbs. weight capacity

> 75AH batteries with 6A charger

> Colours: red, blue or orange











S846 Specifications

Total wt. w/bat. Safe climbo Weight capacity **Battery** Tire Driving dist. Turning radius **Dimensions** LXWXH Motor Speed 56"x 26"x 49" 130kg/286LBS 2.3_{HP} 15крн/9.5мрн 2 х 75ан F12" R12" 50км/32м 161cm/63" 10°/16% 182kg/400lbs All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style