

# Landcruiser 889<sup>XL</sup>SBN **SHOPPRIDER**

## Full-size scooter

Our full-size category provides feature rich luxury in personal mobility. Experience the outdoors in unprecedented style, comfort and performance. This line is for those individuals who want to travel a little further, in style. The full-size series are the most popular scooters in our portfolio.

- > Shoprider 4-wheel scooter
- > Delta tiller with 2-speed (hi-lo)
- > 2HP motor with 12.5 inch tires
- > Full suspension, full light package
- > 75AH batteries with 8A charger
- > Colours: red or silver



**\$5,995.00**  
MSRP



### 889<sup>XL</sup>SBN Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb <sup>o</sup>
227KG/500LBS	56"x 28"x 50"	147KG/323LBS	2.0HP	13.5KPH/8.4MPH	2 x 75AH	4 x 12.5"	48KM/30M	178CM/70"	12°/20%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

[www.eclipsemedical.com](http://www.eclipsemedical.com)

#012023EN